

Product Datasheet

OR11A1 antibody (orb1149769)

Description	Rabbit polyclonal antibody to OR11A1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human OR11A1(NP_039225.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9GZK7
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.