

Product Datasheet

Tau (Phospho-T181) antibody (orb1149825)

Description	Rabbit polyclonal antibody to Tau (Phospho-T181)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T181 of human Tau(NP_001190181.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	74kDa
Uniprot ID	P10636
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.