

Product Datasheet

ARHGEF28 antibody (orb1149847)

Description	Rabbit polyclonal antibody to ARHGEF28
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ARHGEF28 (NP_001073948.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	191kDa
Uniprot ID	Q8N1W1
Entrez	64283
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.