

Product Datasheet

TrkA (Phospho-Y674/675) + TrkB (Phospho-Y706/707) antibody (orb1149848)

Description	Rabbit polyclonal antibody to TrkA (Phospho-Y674/675) + TrkB (Phospho-Y706/707)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y674 & Y675 of human TrkANTRK1(NP_001012331.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	140KDa
Uniprot ID	Q16620 , P04629
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.