

Product Datasheet

TrkB (Phospho-Y516) antibody (orb1149876)

Description	Rabbit polyclonal antibody to TrkB (Phospho-Y516)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y516 of human TrkB (NP_001018074.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q16620
Entrez	4915
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.