

Product Datasheet

GLP1R antibody (orb1149909)

Description	Rabbit polyclonal antibody to GLP1R
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GLP1R (NP_002053.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P43220
Entrez	2740
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.