


Product Datasheet

MLKL (Phospho-S441) antibody (orb1149910)

Description	Rabbit polyclonal antibody to MLKL (Phospho-S441)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S441 of mouse MLKL(NP_001297542.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	54KDa
Uniprot ID	Q9D2Y4
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of C6 ...