

Product Datasheet

Vimentin (Phospho-S83) antibody (orb1149914)

Description	Rabbit polyclonal antibody to Vimentin (Phospho-S83)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S83 of human VIM(NP_003371.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	57KDa
Uniprot ID	P08670
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.