

Product Datasheet

ESR α (Phospho-S167) antibody (orb1149925)

Description	Rabbit polyclonal antibody to ESR α (Phospho-S167)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S167 of human Estrogen Receptor alpha (ESR1)
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	66kDa
Uniprot ID	P03372
Entrez	2099
Dilution Range	WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:100 - 1:200
Expiration Date	12 months from date of receipt.