

Product Datasheet

BEND3 antibody (orb1149953)

Description	Rabbit polyclonal antibody to BEND3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human BEND3 (NP_001073919.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q5T5X7
Entrez	57673
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.