

Product Datasheet

Vinculin antibody (orb1149986)

Description	Rabbit polyclonal antibody to Vinculin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Vinculin (NP_003364.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	124KDa
Uniprot ID	P18206
Entrez	7414
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.