

Product Datasheet

TrkB (Phospho-Y816) antibody (orb1149994)

Description	Rabbit polyclonal antibody to TrkB (Phospho-Y816)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y816 of human TrkB(NP_001018074.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	35kDa/53kDa/59kDa/60kDa/81kDa/91kDa/93kDa
Uniprot ID	Q16620
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.