

Product Datasheet

CMTR1 antibody (orb1150029)

Description	Rabbit polyclonal antibody to CMTR1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CMTR1 (NP_055865.1).
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q8N1G2
Entrez	23070
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.