

Product Datasheet

APP antibody (orb1150153)

Description	Rabbit polyclonal antibody to APP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human APP(NP_000475.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34kDa/72-86kDa
Uniprot ID	P05067
Dilution Range	WB 1:500 - 1:2000IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.