

Product Datasheet

NOP56 antibody (orb1150164)

Description	Rabbit polyclonal antibody to NOP56
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NOP56(NP_006383.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70KDa
Uniprot ID	000567
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.