

Product Datasheet

C19orf48 antibody (orb1150173)

Description	Rabbit polyclonal antibody to C19orf48
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human C19orf48 (NP_954858.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q6RUI8
Entrez	84798
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.