

Product Datasheet

NMDAR1 (Phospho-S897) antibody (orb1150338)

Description	Rabbit polyclonal antibody to NMDAR1 (Phospho-S897)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around of human Phospho-NMDAR1-S897 .
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q05586
Entrez	2902
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.