

Product Datasheet

SLC12A1 antibody (orb1150440)

Description	Rabbit polyclonal antibody to SLC12A1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SLC12A1 (NP_000329.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q13621
Entrez	6557
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.