

Product Datasheet

SARS-CoV-2 ORF6 antibody (orb1150455)

Description	Rabbit polyclonal antibody to SARS-CoV-2 ORF6
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of SARS-CoV-2 ORF6.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	12KDa
Uniprot ID	P0DTC6
Entrez	43740572
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.