

Product Datasheet

TRAPPC1 antibody (orb1150566)

Description	Rabbit polyclonal antibody to TRAPPC1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human TRAPPC1(NP_067033.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5R8
Dilution Range	IHC 1:50 - 1:100
Expiration Date	12 months from date of receipt.