

Product Datasheet

BrdU antibody (orb1150570)

Description	Rabbit polyclonal antibody to BrdU
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	BrdU
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC 1:50 - 1:500 IF 1:50 - 1:500
Expiration Date	12 months from date of receipt.