

Product Datasheet

SMAD2/SMAD3 antibody (orb1150581)

Description	Rabbit polyclonal antibody to SMAD2/SMAD3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated of human SMAD2/SMAD3SMAD3(NP_005893.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52kDa/60kDa
Uniprot ID	P84022, Q15796
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.