

Product Datasheet

SARS-CoV-2 NSP1 antibody (orb1150608)

Description	Rabbit polyclonal antibody to SARS-CoV-2 NSP1
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of coronavirus NSP1(YP_009742608.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20KDa
Uniprot ID	P0DTC2
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.