

Product Datasheet

GSTA4 antibody (orb1150705)

Description	Rabbit polyclonal antibody to GSTA4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GSTA4 (NP_001503.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	26kDa
Uniprot ID	O15217
Entrez	2941
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.