

Product Datasheet

Acetyl-Histone H4-K5/K8/K12/K16 antibody (orb1150726)

Description	Rabbit polyclonal antibody to Acetyl-Histone H4-K5/K8/K12/K16
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic acetylated peptide around K5/K8/K12/K16 of human HistoneH4.
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	11kDa
Uniprot ID	P62805
Dilution Range	IHC 1:50 - 1:200
Expiration Date	12 months from date of receipt.