

Product Datasheet

PIEZO1 antibody (orb1150766)

Description	Rabbit polyclonal antibody to PIEZO1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human FAM38A/PIEZO1(NP_001136336.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	286kDa
Uniprot ID	Q92508
Dilution Range	WB 1:500 - 1:2000IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.