

Product Datasheet

Tau (Phospho-T548) antibody (orb1150800)

Description	Rabbit polyclonal antibody to Tau (Phospho-T548)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around of human Phospho-Tau-T548.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P10636
Entrez	4137
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:200
Expiration Date	12 months from date of receipt.