

Product Datasheet

TrkB (Phospho-Y705) antibody (orb1150807)

Description	Rabbit polyclonal antibody to TrkB (Phospho-Y705)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y705 of human Phospho-TrkB-Y705
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	140kDa
Uniprot ID	Q16620
Entrez	4915
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.