

Product Datasheet

SARS-CoV-2 3CLpro antibody (orb1150829)

Description	Rabbit polyclonal antibody to SARS-CoV-2 3CLpro
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of coronavirus 3CLpro(YP_009725295.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30-35KDa
Uniprot ID	P0DTC2
Dilution Range	WB 1:500 - 1:2000 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.