

Product Datasheet

STING/TMEM173 antibody (orb1150899)

Description	Rabbit polyclonal antibody to STING/TMEM173
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human STING/TMEM173 (NP_938023.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	35-40KDa
Uniprot ID	Q86WV6
Entrez	340061
Dilution Range	WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.