

Product Datasheet

TBC1D5 antibody (orb1150967)

Description	Rabbit polyclonal antibody to TBC1D5
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TBC1D5(NP_055559.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	91kDa
Uniprot ID	Q92609
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.