

Product Datasheet

SARS-CoV-2 Nucleoprotein antibody (orb1151059)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Nucleoprotein
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	DOT, ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of coronavirus Nucleoprotein(P0DTC9).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Dilution Range	WB 1:500 - 1:2000IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.