

Product Datasheet

DDR1 (Phospho-Y513) antibody (orb1151071)

Description	Rabbit polyclonal antibody to DDR1 (Phospho-Y513)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y513 of human DDR1 (NP_054700.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q08345
Entrez	780
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.