

Product Datasheet

ASB17 antibody (orb1151072)

Description	Rabbit polyclonal antibody to ASB17
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of mouse ASB17(NP_080034.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34KDa
Uniprot ID	Q8VHP9
Dilution Range	WB 1:1000 - 1:2000
Expiration Date	12 months from date of receipt.