

Product Datasheet

ALDH2 antibody (orb1151082)

Description	Rabbit polyclonal antibody to ALDH2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human ALDH2(NP_000681.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50-56KDa
Uniprot ID	P05091
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.