

Product Datasheet

SURF4 antibody (orb1151084)

Description	Rabbit polyclonal antibody to SURF4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide of human SURF4 (NP_001267719.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	30KDa
Uniprot ID	015260
Entrez	6836
Dilution Range	WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.