

Product Datasheet

TRIM10 antibody (orb1151087)

Description	Rabbit polyclonal antibody to TRIM10
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TRIM10 (NP_439893.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9UDY6
Entrez	10107
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.