

Product Datasheet

IL18 antibody (orb1151088)

Description	Rabbit polyclonal antibody to IL18
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human IL18(NP_001553.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18KDa/22KDa
Uniprot ID	Q14116
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.