

Product Datasheet

SLC22A11 antibody (orb1151089)

Description	Rabbit polyclonal antibody to SLC22A11
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SLC22A11(NP_060954.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65KDa
Uniprot ID	Q9NSA0
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.