

Product Datasheet

UGT2B17 antibody (orb1151090)

Description	Rabbit polyclonal antibody to UGT2B17
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human UGT2B17(NP_001068.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61kDa
Uniprot ID	075795
Dilution Range	IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.