

Product Datasheet

Pan-Crotonyl-K antibody (orb1151104)

Description	Rabbit polyclonal antibody to Pan-Crotonyl-K
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to a sequence containing acetylatedK.
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	10-250KDa
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.