

Product Datasheet

MEIS3 antibody (orb1151106)

Description	Rabbit polyclonal antibody to MEIS3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MEIS3(NP_001287988.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	Q99687
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.