

Product Datasheet

MSY2 antibody (orb1151107)

Description	Rabbit polyclonal antibody to MSY2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human MSY2(NP_057066.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50KDa
Uniprot ID	Q9Y2T7
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.