

Product Datasheet

HASPIN antibody (orb1151110)

Description	Rabbit polyclonal antibody to HASPIN
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HASPIN(NP_114171.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100KDa
Uniprot ID	Q8TF76
Dilution Range	WB 1:1000 - 1:2000
Expiration Date	12 months from date of receipt.