

Product Datasheet

Methyl-HDAC3-H134/H135 antibody (orb1151113)

Description	Rabbit polyclonal antibody to Methyl-HDAC3-H134/H135
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human Methyl-HDAC3-H134/H135, 3-methyl (tau-methyl).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O15379
Entrez	8841
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.