

Product Datasheet

ZNF75A antibody (orb1151114)

Description	Rabbit polyclonal antibody to ZNF75A
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human ZNF75A(NP_694573.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35kDa
Uniprot ID	Q96N20
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.