

Product Datasheet

2'-O-Methylguanosine(Gm) antibody (orb1151115)

Description	Rabbit polyclonal antibody to 2'-O-Methylguanosine(Gm)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Immunogen	Chemical compounds corresponding to 2'-O-Methylguanosine(Gm).
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.