

## Product Datasheet

### Histone H2A.1 antibody (orb1151116)

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Histone H2A.1                     |
| <b>Reactivity</b>          | All, Human, Mouse, Rat                                          |
| <b>Conjugation</b>         | Unconjugated                                                    |
| <b>Tested Applications</b> | IHC, WB                                                         |
| <b>Immunogen</b>           | Recombinant fusion protein of human Histone H2A.1(NP_064707.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal,50% glycerol,pH7.3.                   |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                     |
| <b>Note</b>                | For research use only                                           |
| <b>Isotype</b>             | IgG                                                             |
| <b>Clonality</b>           | Polyclonal                                                      |
| <b>MW</b>                  | 14kDa                                                           |
| <b>Uniprot ID</b>          | <b>P0C0S8</b>                                                   |
| <b>Dilution Range</b>      | WB 1:500 - 1:2000IHC 1:50 - 1:200                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                 |