

Product Datasheet

GLYATL1 antibody (orb1151122)

Description	Rabbit polyclonal antibody to GLYATL1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human GLYATL1 (NP_001207423.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	35kDa
Uniprot ID	Q969I3
Entrez	92292
Dilution Range	WB 1:500 - 1:2000 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.