

Product Datasheet

C1QTNF3 antibody (orb1151123)

Description	Rabbit polyclonal antibody to C1QTNF3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human C1QTNF3 (NP_112207.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9BXJ4
Entrez	114899
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.